

LabSpec 5 (LS 5)



LabSpec 5



Config LabSpec
5

HORIBA JOBIN YVON

Explore the future

HORIBA

Notation

- Names of icons and menu items are denoted in **bold**.
- Names of windows are denoted in *italic*.
- Names of columns, rows or entries in a window are denoted in all CAPITAL.

Configuration Utility

HORIBA JOBIN YVON

Explore the future


HORIBA

Versatility


- We provide a wide range of Raman microscope/spectrometer instruments (analytical, research, processing, etc.), which require various configurations to meet the need of the applications specified by customers.
- While we can optimize the system to yield the best results in the target application, the service (maintaining existing systems, adding accessories, upgrading obsolete parts, etc.) could become a daunting task.
- The design goal of LabSpec 5 is to provide a platform that can adapt and control dynamically changing configurations without having to replace the software.
- The task is accomplished mainly through Config utility, where the specific configuration of the instrument is specified.

Accessing Config

- Double clicking **Config** icon opens the *Config Login* window, the operator can access the configuration of the LabSpec 5 in three levels.
 - No password: Clicking **OK** button with no password will prompt a message box that informs that changing configuration will be disabled. Clicking **OK** button of the message box opens the *Config* window in read-only format.



Config



- Each element comes with a list menu, from which the specific component and its configurations can be defined.
- Module option is where customized options are listed. When activated, corresponding icons or controls will become available in LabSpec 5 main window.

Overview

HORIBA JOBIN YVON

Explore the future

HORIBA

Main Window

The screenshot shows the LabSpec software interface. At the top is a menu bar with 'File', 'Edit', 'Data', 'Options', 'Acquisition', 'Setup', 'Maintenance', 'Scripts', and 'Help'. Below the menu bar is a toolbar with various icons. The main workspace is a large grey area. At the bottom is a control panel with several sections: 'Laser' (with a dropdown menu showing 'External La'), 'Filter' (with a dropdown menu), 'Hole' (with a value of '115 μm'), 'Slit' (with a value of '117 μm'), 'Spectrometer' (with a value of '277978 cm⁻¹'), 'Acquisition' (with three input fields containing '1', '1', and '1'), and 'Aramis' (with a 'Setup' button). A 'STOP' button is located in the top right corner of the window.

Callouts in the image point to the following features:

- File I/O
- Experiment setup
- Data acquisition
- Axes and Cursors
- Data Information
- Data processing
- Switch lasers
- Neutral density filter
- Confocal hole diameter
- Slit width
- Spectrometer
- Grating/Objective
- Exposure time
- Laser block/unblock

HORIBA For Help, press F1

Video Acquisition

The screenshot displays the LabSpec software interface. At the top, the menu bar includes File, Edit, Data, Options, Acquisition, Setup, Maintenance, Scripts, and Help. The toolbar contains various icons, with the Video Acquisition icon (a camera) circled in red. A white text box with a black border is positioned over the main workspace, containing the text: "Video Continuously capture and display the video image Both camera and lamp can be polarized." An arrow points from the circled icon to this text box. Another arrow points from the text box to the Polarization dialog box. The Polarization dialog box is open, showing settings for the White Lamp and Camera. Each section has two buttons for "0 deg" and "90 deg", and a text input field with "179" and a label "0 -> 180 deg". At the bottom of the interface, there are several control panels: Laser (External La), Filter, Hole (115 μm), Slit (117 μm), Spectrometer (277978 cm⁻¹), Acquisition (1, 1, 1), and Aramis (Setup). A footer note says "For Help, press F1".

Video

The screenshot displays the LabSpec software interface. The main window is titled "LabSpec - Video : ch-5 map2". The menu bar includes File, Edit, Data, Options, Acquisition, Setup, Maintenance, Scripts, Window, and Help. A toolbar with various icons is located below the menu bar. On the left side, a vertical toolbar contains icons for navigation and zooming. A white text box with a bracket points to this toolbar, containing the text "Menu options available for Video".

The central area features a video feed window titled "Video : ch-5 map2". The video shows a grayscale image of a sample with a coordinate system. The X-axis is labeled "X (μm)" and ranges from -100 to 100. The Y-axis is labeled "Y (μm)" and ranges from -100 to 100. A green scale bar in the bottom right corner of the video indicates "20 μm".

At the bottom of the interface is a control panel with several sections:

- Laser:** External La
- Filter:** [Dropdown menu]
- Hole:** 115 μm
- Slit:** 117 μm
- Spectrometer:** 277978 cm⁻¹
- Acquisition:** 1800, x100, Test
- Aramis:** [Buttons for Start, Stop, Setup]

At the bottom left, there is a text prompt: "For Help, press F1". At the bottom right, there are status indicators: "X: 174.084 Y: 129.891 I: 57".

HORIBA

Explore the future

HORIBA

Spectrum RTD (Real Time Display)


The screenshot displays the LabSpec software interface. At the top, the menu bar includes File, Edit, Data, Options, Acquisition, Setup, Maintenance, Scripts, and Help. A toolbar below the menu bar contains various icons, with the Spectrum RTD icon (a circular icon with a spectrum) circled in red. A callout box points to this icon, stating: "Spectrum RTD Continuously capture and display the spectrum at the current spectrometer position".

Another callout box points to the Spectrometer section of the bottom control panel, stating: "Spectrometer position The spectral range is defined by the size of the CCD, laser wavelength and grating".

A third callout box points to the Acquisition section of the bottom control panel, stating: "RTD exposure time is set separately from the Data Acquisition exposure time".

The bottom control panel includes sections for Laser (External La), Filter, Hole (115 μm), Slit (117 μm), Spectrometer (277978 cm⁻¹), Acquisition (1, 1, 1), and Aramis (Setup). A status bar at the bottom left reads "For Help, press F1".

Spectrum



Spectrum Acquisition

The screenshot displays the LabSpec software interface. The menu bar includes File, Edit, Data, Options, Acquisition, Setup, Maintenance, Scripts, and Help. The toolbar contains various icons, with two icons circled in red: a multi-range icon and a single range icon. Arrows point from these icons to two text boxes:

- Extended Range Acquisition**
Setup the spectral range(s) for the **Spectrum Acquisition**
- Spectrum Acquisition**
Acquire a single spectrum

The 'Extended range' dialog box is open, showing the following settings:

From	To	Time(%)	
300	380	100	<input type="checkbox"/>
400	480	100	<input type="checkbox"/>
500	580	100	<input type="checkbox"/>
600	680	100	<input type="checkbox"/>
700	780	380	<input type="checkbox"/>

Additional settings in the dialog box:

- Min overlap(pix): 100
- SubPixel: 1
- Use extended range:
- Auto overlap:
- Combine data:
- Adjust intensity:
- Return to first window:


Buttons at the bottom of the dialog box: OK, Cancel, Help.

The bottom status bar of the LabSpec window shows the following parameters:

- Laser: External La
- Filter: [Icon]
- Hole: 115 μm
- Slit: 117 μm
- Spectrometer: 277978 cm^{-1}
- Acquisition: 1800, x100, Test
- Araxis: [Icon]

For Help, press F1

Spectrum



Map Acquisition

The screenshot shows the LabSpec software interface for map acquisition. The main window is titled "LabSpec - Video : Crack1" and contains a menu bar (File, Edit, Data, Options, Acquisition, Setup, Maintenance, Scripts, Window, Help) and a toolbar with various icons. A vertical toolbar on the left side is annotated with a white box containing the text: "Interactive menu options to setup (lateral) map." The central area displays a video feed of a crack, with a coordinate system overlaid. The X-axis is labeled "X (μm)" and ranges from -50 to 50. The Y-axis is labeled "Y (μm)" and ranges from -40 to 40. A scale bar in the bottom right of the video window indicates "10 μm". The bottom of the interface features a control panel with several sections: "Laser" (External La), "Filter", "Hole" (115 μm), "Slit" (117 μm), "Spectrometer" (277978 cm⁻¹), "Acquisition" (1800, x100, Test), and "Aramis" (Setup). The status bar at the bottom shows "For Help, press F1" and numerical values "46.03" and "34.95".

HORIBA LABSPEX

Explore the future

HORIBA

Map Acquisition



Mapping Properties
Setup the map and/or batch acquisition

The Mapping Properties dialog box contains a table with the following data:

	Array Size	From	To	Increment	Formula	Unit	elati
Time	<input type="checkbox"/>	11	0	100.0	1.000	<input type="checkbox"/>	sec
Custom	<input type="checkbox"/>	11	0	100.0	1.000	<input type="checkbox"/>	
Y	<input type="checkbox"/>	5	-1293	1293	1.000	<input type="checkbox"/>	
X	<input type="checkbox"/>	10	-1701	1701	1.000	<input type="checkbox"/>	

Below the table is a checkbox labeled "Always show XY parameters" which is currently unchecked. On the right side of the dialog are buttons for "OK", "Apply", "Help", and "Cancel".

The Devices dialog box shows a list of device parameters:

- Custom
- Time
- Exposure time
- Accumulation
- Gain
- Slit
- Spectro
- Hole
- Premono2
- Shutter1
- Shutter2
- Line/Point
- Video

Buttons for "OK" and "Cancel" are on the right.

The hardware configuration panel includes the following sections:

- Laser:** External La
- Filter:** [Color filter icon]
- Hole:** 115 μm
- Slit:** 117 μm
- Spectrometer:** 277978 cm^{-1}
- Acquisition:** 1800, x100, Test
- Aramis:** Setup

For Help, press F1

Map Acquisition

LabSpec

File Edit Data Options Acquisition Setup Maintenance Scripts Help

Mapping Acquisition

Acquire a series of spectra as setup in **Mapping Properties** and **Extended Range Acquisition**

External La

Filter

Hole 115 μm

Slit 117 μm

Spectrometer 277978 cm^{-1}

1800

x100

Test

Acquisition

1

1


1

Aramis

Setup

For Help, press F1

Map



Max/min spectral intensity across the entire data set

Video image overlaid with the integrated intensity maps

Spectrum at the current cursor position in the **Video** or **Map**

Integrated intensity maps defined by three pairs of cursors (color coded)

Maintenance

HORIBA JOBIN YVON

Explore the future



HORIBA

Spectral Axis Calibration

- X-axis (spectral axis) of a Raman microscope is calculated using the full trigonometric relationship involving the included spectrometer angle.
- When the instrument is well aligned, there is a linear relationship between the spectrometer parameters and the spectral axis coordinates.
- The slope (KOEFF) and intercept (ZERO) of this linear relationship can be 'tweaked' to adapt to the daily fluctuation of the environment (e.g. temperature).
- The best practice is to keep record of changes made to these parameters.
- Please contact HORIBA Jobin Yvon if a significant drift in spectral axis is suspected or observed.

Spectral Axis Calibration

1. Place Si reference target on the sample stage
2. Change the X-axis unit to nm by selecting **Options / Unit / nm**
3. Set Spectrometer to **Zero Order**



4. Check zero position (laser line) using **Spectral Real Time Display**
5. If the laser line position is off zero, select **Setup / Instrument Calibration** and adjust the value for ZERO. Make sure to check zero position after every adjustments by repeating steps 3 and 4.
6. Change the X-axis unit back to 1/cm by selecting **Options / Unit / 1/cm**
7. Set Spectrometer to 520 cm⁻¹, where Si line should appear.
8. Check Si line position using **Spectral Real Time Display**
9. If the laser line position is off 520, select **Setup / Instrument Calibration** and adjust the value for KOEFF. Make sure to check line position after every adjustments by repeating steps 7 and 8.

Spectral Axis Calibration

- To ensure the linearity over the full spectral range, the procedure in the previous page can be repeated with different spectrometer position.
- From the step 7, set the spectrometer to 600, 700, etc., and repeat the rest of the steps.
- Sharp and strong peaks from other samples of well known spectra can be used as well. (e.g. Aspirin, Chloroform)


Data Acquisition

HORIBA JOBIN YVON


Explore the future

HORIBA


Video Acquisition – Real Time



Video
Continuously
capture and display
the video image



Both light source
(white lamp)
and detector (camera)
can be polarized



Video

The screenshot displays the LabSpec software interface. The main window is titled "LabSpec - Video : ch-5 map2" and contains a menu bar (File, Edit, Data, Options, Acquisition, Setup, Maintenance, Scripts, Window, Help) and a toolbar. A vertical toolbar on the left side is highlighted with a white box and labeled "Menu options available for Video".

In the center, a video window titled "Video : ch-5 map2" shows a grayscale image of a sample. The image has a coordinate system with X and Y axes in micrometers (μm), ranging from -100 to 100. A green scale bar in the bottom right corner indicates 20 μm .

At the bottom, there are several control panels:

- Laser:** External La
- Filter:** [Dropdown menu]
- Hole:** 115 μm
- Slit:** 117 μm
- Spectrometer:** 277978 cm^{-1}
- Acquisition:** 1800, x100, Test
- Aramis:** Setup

At the bottom left, it says "For Help, press F1". At the bottom right, it shows coordinates: X: 174.084, Y: 129.891, I: 57.


HORIBA

Explore the future


HORIBA

Transmission Light Pyroxine

Parallel polarization




Cross polarization




Reflection Light



Parallel polarization



Cross polarization




Extended Video Acquisition




- Requires a motorized xy-stage.
- Select NGSExtendedVideo.dll in *Config* to activate, which adds an additional icon to the LS 5 tool bar
- Selecting the icon opens the *Extended Video* window for setting up the extended video acquisition.

HORIBA JOBIN YVON


Extended Images



- Extending the field of view by patching multiple images together.
- This is **NOT** a real time image.



Move to the Cursor Position




- Requires a motorized xy-stage.
- Select NGSMouseTable.dll in *Config* to activate, which adds additional icons to the left tool bar of **Video**
- When the **Target** is selected, a bulls-eye appears on the video image, which defines the current center.
- Clicking on a new location in the image moves the bulls-eye and the stage moves automatically to center at the new position.

Spectrum Acquisition – Real Time



Spectrum RTD ← → Detector Image RTD



- **Acquisition / RTD**

- There are two Real Time Display (RTD) modes
- Selecting **Acquisition / RTD** menu opens RTD window
- Image binning X and Image binning Y apply to **Detector Image RTD**.
 - Number of pixels specified here are added and read out as one pixel.
 - Increasing number of binning speeds up the read-out process.
- Readout zone applies to both **Detector Image RTD** and **Spectrum RTD**.
 - There are three Readout zone modes: **Spectrum**, **Image** and **Full**
 - **Spectrum** and **Image** are defined in **Acquisition / Detector** and window the detector read-out regions
 - Usually, Readout zone is set to **Spectrum**
- When selected, Accumulate image takes the number of accumulations (set in the bottom tool bar) into account.
- (RTD) Exposure time can be modified here as well as in the bottom tool bar.

- RTD is triggered on by clicking Icons on the LS 5 tool bar.


Spectrum Acquisition - Options




Spectrum RTD ← → Detector Image RTD

- **Acquisition / Options**


- Exposure time and Accumulation number can be set here as well as in the bottom tool bar
- The name of the spectrum (or the data set) is set to Data name and a increasing integer suffix by default
- Refresh time: Refresh the data measurement in progress at the set interval
- Binning factor: Read out the set number of pixels as one. Increasing binning factor lowers the spectral resolution while increasing the signal.
- Accumulation mode: Average, Sum, Detector
- Spike filter: None, Multi (auto add), Multi, Single
- Autofocus: Off, Before acquisition, Before each point
- Scanning: Stage (point mapping mode), Scanner (line scanning mode).
- Use detector : Active detector, Both
- Signal mode : Signal, Signal-Dark




Extended Range Acquisition



- Selecting **Acquisition / Extended Range** menu or **Extended Range** icon opens *Extended Range* window
- Select Auto-overlap and Combine Data options to use CREST (Continuous Rapid Extended Scanning Technology) to measure the spectrum.
- Spectrum acquisition is started by clicking **Spectrum Acquisition**.



Map Acquisition



LabSpec - Video : Crack1

File Edit Data Options Acquisition Setup Maintenance Scripts Window Help

Mapping Acquisition ← → Mapping Properties

- Setting up the surface map
 - Icons to define mapping properties appear on the left tool bar when the *Video* window is selected.
 - Each icon represents a map of a different shape.
 - Selecting an icon overlays the selected shape on top of the video image (pink dots).
 - The size of the map is adjusted by clicking and dragging with the mouse.


Map Acquisition



Mapping Acquisition



Mapping Properties

- Setting up the surface map
 - **Mapping properties** icon opens *Mapping properties* window
 - Array Size and Increment must be defined here.




Map Acquisition – Other Parameters

- *Mapping properties* window provide the place to set up parameters other than x and/or y scanning.
 - Depth profile is performed by employing the parameter Z
 - Other variables such as slit width and confocal hole diameters can be set to scan over a range at an increment.
 - Right clicking on the parameter list activates context menu. Select **Insert** to open *Device* window. Select variables to scan. (e.g. SLIT, HOLE)



Map Acquisition – Other Parameters

- Priorities of the variables
 - The order of the variables sets their priority during the mapping. The bottom parameter has higher priority, the top parameter has lower priority.
 - When set as below, the measurements will be made along x-axis at a constant y position before moving to a new y position to repeat the process.
 - To change the order of variables, right click on the variable name to activate the context menu. Use **Remove**, **Up**, **Down** menus to change the list and order of variables.




HORIBA JOURNAL VOIN

Map Acquisition – Options

- Increment
 - When the INCREMENT box is checked, the measurements are performed at the set increment.
 - When the INCREMENT box is not checked, measurements are performed at the increment calculated from FROM, TO and ARRAY SIZE.
- Relative
 - When the RELATIVE box is checked, then the current position of the image is the (0,0,0) position of the map table.
 - When the RELATIVE box is not checked, then the (0,0,0) position of the image is the (0,0,0) position of the map table.

Map Acquisition – Options

- Extra Images
 - **Acquisition / Extra images**
 - When the autofocus is activated, extra images like the Reflection (on the autofocus detector) and the Z position (after each autofocusing) are available.
 - These are recorded
 - At the same time as the Raman mapping if you select MAP
 - Without recording any spectrum if you select IMAGE and click on IMAGE.
 - If you use IMAGE then you only obtain a topographic picture of the sample and the variation of reflectivity over the surface of the sample.



Data Treatment


HORIBA JOBIN YVON

Explore the future

HORIBA

Basic Data Treatment

- Baseline correction
- Spectral normalization and correction
- Smoothing
- Fourier transform
- Math
- Peak searching and fitting
- Profile from image







Baseline Correction

Build the baseline by adding points at the cursor position

Remove the points defining the baseline



Build the baseline automatically and subtract from the current spectrum


Build the baseline automatically but NO subtraction

Copy the baseline

Paste the baseline to an other spectrum




Spectral Normalization and Correction



- **Normalize:** Normalizes all traces to same area value (100).
- **Zero:** Moves all traces to the minimum intensity level.
- **Get:** Takes the activated spectrum and put it in Corrector frame. This data object will be used as parameter.
- **Delete:** Removes the Corrector spectrum.
- **Subtract:** Subtracts the Corrector spectrum from all traces.
- **Multiply:** Multiplies the trace on the Corrector spectrum.
- **Correct:** Subtracts the Corrector spectrum from all traces. The intensity of the Corrector is multiplied to fit the intensity of the trace
- **Thr:** Remove all trace with maximum intensity less then threshold value (in % to maximum data intensity)



Spectral Normalization and Correction




LIMITS: Check the box if you want that only a part of the trace is used for operation. This parameter is used for **Normalize** and **Correct** operation.

The close-up shows the 'Limits' section of the software interface. It contains two input fields: 'From' with the value 506.372 and 'To' with the value 1551.65. To the right of these fields is a checked checkbox. Red arrows point from the text above to the 'From' field, the 'To' field, and the checkbox.

FROM and TO: Edit these fields to modify the limits of the selected region (linked to the red vertical cursors)



Smoothing



Savitsky-Golay smoothing.


First derivative

Second derivative

Median filter smoothing




Fourier Transform



- This is another smoothing tool
- When a spectrum is Fourier transformed, the low frequency data (spectral features) and high frequency data (noise) are separated.
- Therefore, eliminating high frequency region in the Fourier transform domain removes high frequency noise.
 - Fourier transform the spectrum
 - Cut off high frequency region
 - Reverse Fourier transform to spectrum



Fourier Transform




None No apodization function
Line Linear function
Sqr Parabolic function
Cos Cosine function

None: No filter function
Traffic: Traffic function

- LIMIT: Position of the cut off point in %. This parameter defines the smoothing factor.
 - 0: Fully smoothed,
 - 100: No smoothing.
- Click on the **Go** button to apply the value or drag the pointer to select the LIMIT visually.
- APOD: Type of the apodization function with apodization value reaching zero at a Limit point.
- FILTER: Select here the type of the filter function



Math



Combine data objects into a single object

- FUNC 1 and FUNC 2: Basic math functions (+, -, *, /, ^), negating, power, log, exp, sin, cos, asin, acos, atan, abs, sqrt, step.
- **Go**: Perform the operation
- x : denotes the active object
- y : denotes other objects of the active window
- a, b, c, d..... : values of the axes for each data points (next slide)



Examples

If we apply the formula ' $x+2*a$ ' to a spectrum, x is the intensity and a the Raman at each point.

Before operation

Raman shift	100	200	300	...	a
Intensity	200	250	500	...	x

After operation

Intensity values	20000	50000	150000	...
------------------	-------	-------	--------	-----

If we apply the formula ' $x+b-a$ ' to an image, x is the intensity, a is the x distance and b is the y distance

Before operation



	10	20	30	...	a
40	500	250	500		
50	50	120	400		
60	30	0	500		
...					
b					

After operation

	530	270	510
	90	150	420
	80	40	530



Math




Operations using constants can be performed interactively using icons as well.

Add/subtract a constant

Multiply/divide by a constant




Peak Searching and Fitting



Insert a peak

Adjust position,
height and width

Remove the selected
peak



Search: search peaks automatically

Approx: approximate the peak position, height and width without fitting

Fit: start iteration

Init: initialize

Clear: clear all peaks and calculation results if there is any

Convert: convert the fitted spectrum with the calculation result

Peaks: open *Peak* window

Variables: open *Variables* window



Peak Searching and Fitting

Peaks

	p	Min	Max	Fix	a	Fix	w	Fix	g	Fix	Formula
1	512.948	511.26	515.26	<input type="checkbox"/>	102.373	<input type="checkbox"/>	16.2417	<input type="checkbox"/>	1	<input type="checkbox"/>	GaussLoren()
2	524.881	523.73	527.73	<input type="checkbox"/>	38.2526	<input type="checkbox"/>	20.8675	<input type="checkbox"/>	1	<input type="checkbox"/>	GaussLoren()
3	616.472	616.472	620.472	<input type="checkbox"/>	0	<input type="checkbox"/>	0.001	<input type="checkbox"/>	1	<input type="checkbox"/>	GaussLoren()
4	636.465	633.617	637.617	<input type="checkbox"/>	204.964	<input type="checkbox"/>	26.9624	<input type="checkbox"/>	0	<input type="checkbox"/>	GaussLoren()

Remove row
Insert row

OK
Copy
Paste
Apply
Close

Right clicking in the table activate the context menu, which allows add or remove peaks to fit.


Parameters of fitting functions can be adjusted for each peak from *Variables* window

The entire peak fitting table can be copied and pasted using the **Copy** and **Paste** button.


Variables

	Init	Min	Max	Copy
p	x	p-2	p+2	<input checked="" type="checkbox"/>
w	1	0.001	100	<input type="checkbox"/>
a	y	0	maxy	<input type="checkbox"/>
b	0	-maxy	maxy	<input type="checkbox"/>
c	0	miny	maxy	<input type="checkbox"/>
d	0	-maxy	maxy	<input type="checkbox"/>
g	0.5	0	1	<input type="checkbox"/>

OK
Cancel



Profile from Image




Vertical or horizontal profile over an image can be obtained from Profile window.

You can apply to a profile the same operations as for a spectrum, then using the **Corr** button to put the corrected data back to the image where the profile was extracted

Advanced Data Treatment

- Analyzing map image and profiles
- New spectral profile

Analyzing Images and Profiles



HORIBA
For Help, press F1

Explore the future

HORIBA



Analyzing with Cursors

Check USE box to generate the corresponding map.

Check BASE box to subtract baseline.

Check GREEN/BLUE box to generate the map of the corresponding ratio.

Check SPECTRUM box to see the spectrum of the selected data point(s)

Options	Use	Base	From	To
Red	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	919.462	2147.34
Green	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1297.04	1355.36
Blue	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	465.146	560.307

Green/Blue

Spectrum

Correct

OK Convert Help Cancel

Edit FROM and TO fields to adjust the limits of the spectral regions.

Click on the **Correct** button to modify data on the selected data point(s).



Analyzing with Cursors


The initial method for analysing a mapped image or profile is to use three cursors (**R**, **G** and **B**) to define regions.

The image is generated by displaying the spectral intensity between these cursors.

The three cursors are simply selected from the three coloured cursor.

These three cursors appear when the window Splm is activated.


If one or more of the cursors aren't currently displayed on screen, they can be returned by clicking on the **Cursor normalization** icon on the toolbar.



Example

- Sample description
 - A raw pharmaceutical pellet. The map shows the distribution of 3 widespread fillers in pharmaceutical pellets: Cellulose, Starch and Magnesium stearate.
- Measurement
 - A large map is recorded onto the surface of the pellet.
- Acquisition conditions
 - Instrument: LabRam HR
 - Laser wavelength : 633 nm
 - Grating : 600 gr/mm
 - Confocal hole : 600 μm
 - Objective: $\times 50$
 - Map size: 400 \times 400 μm and 36 \times 36 data points
 - Acquisition time: 10 s/point

The map displays the distribution of 3 components: Starch, Cellulose, and MgStearate




HORIBA JOBIN YVON

Explore the future

HORIBA

To display the integrated intensity map, position pairs of Red, Green, and Blue cursors around selected Raman bands in *Point* or *SpIm* by selecting the corresponding color arrow on the left toolbar.




HORIBA JOBIN YVON

Explore the future

HORIBA

Sometimes, an intensity profile along a line can provide useful information such as diameter of agglomerations. Select a map (right tool bar) and click **Profile** icon (top tool bar). The intensity profile of the selected intensity map is displayed along the vertical or horizontal cursor line (red cross in *Map* window).



Select intensity map

HORIBAJOBIN YVON

Explore the future

HORIBA



Modeling

Model

Properties

Name

Factors

Thr(%)

Allow negative scores

Show error map

Show model sum

Normalize model

Get

Create

Help

Close

Supervised or unsupervised factors

Normalize model

If the spectra of the pure components are available, they can be taken into the modeling program.

If not, LabSpec 5 allows you to create a model automatically by using a factor analysis algorithm.

In terms of normalization of the reference spectra:


If the original spectra had not been recorded in the same conditions you wouldn't lose any information by normalizing because the relative intensities were already unreliable.

But, if the spectra have been recorded in the same conditions, the relative intensities are representative of the relative Raman scattering power of each component. In this latter case, you lose information by normalizing.


Example

- Sample description
 - AsO_4^{3-} can substitute for SO_4^{2-} in sulfate minerals. Raman spectrometry proved very useful for investigating the formation of SO_4^{2-} - AsO_4^{3-} intermediate phases.
 - The sample comes from a former industrial site in France comprising mainly sulfate minerals similar to gypsum ($\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$) and thenardite (Na_2SO_4) and with a high As content (8 to 18% w/w).
- Acquisition conditions
 - Instrument: LabRam HR
 - Laser wavelength : 633 nm
 - Map size: 51 × 70 μm
 - Step size: 1 μm
 - Acquisition time: 10 s/point

Two phases are observed. The characteristic peak of each phase is displayed between Red and Green cursors:



To extract two pure spectra constituting the sample, click **Modeling** icon to open *Model* window. Set FACTORS to 2 and click **Create**.



HORIBA

For Help, press F1

S: 988.043 S: 1026.3 P: 38.2531


Start | Eunah Lee - Inbox - Lotu... | Microsoft PowerPoint - [L... | LabSpec - SpIm : rosa...

98% | 8:59 PM


Explore the future

HORIBA


Two phases are extracted by the Modeling procedure. The characteristic peak of each pure phase is displayed in Red and Green.




MODE: Overlay






To change the display mode of score images, right click in the image to activate window context menu and select **Format and scale**. Select **TILE** from MODE.




MODE: Tile




Open the corresponding video image, right click in the image to activate the window context menu, and select **Imposition**. Check OVERLAY IMAGES box. Select an intensity map or a score map to overlay.





It is also possible to use the band fitting results to create a map of peak position (p), amplitude (a), band width (w) and the area (a) of the peak.

Check MAP box of the parameter to map and click **Apply**.




Example

- Sample description
 - 1 mm square silicon structure.
- Measurements
 - A line measurement is recorded at mid-height of the square from one edge to the other.
 - A spectrum is recorded at every 15 μm .
 - The Silicon peak position in a non-stressed Si crystal is 520.7cm^{-1} .
 - If a tensile or a compressive stress is occurring in the structure (due to a defect, a coating, a boundary, etc.), then the position of the peak shifts to lower (tensile) or higher (compressive) wavenumbers.
- Acquisition conditions
 - Instrument: LabRam HR
 - Laser wavelength : 532 nm
 - Grating: 1800 gr/mm
 - Confocal hole diameter: 200 μm
 - Step size (x-axis only): 15 μm
 - Extended video image: the picture is recorded over a very large surface area.

SpIm displays all the spectra recorded along the line. *Point* shows the spectrum at the cursor position selected either in *Video* or *Map*. The line of analysis is superimposed on the video image. *Map* displays up to three parameters that are selected in *SpIm*.

Red cursors are position in *SpIm* to map the intensity of the silicon peak along the line. The intensity drops at the edge of the silicon square.





HORIBA



Explore the future

HORIBA


To display *SpIm* in 2D, right click in *SpIm* window to activate the window context menu and select **Format and scale**. Select 2D from DIMENSION



To display the line of measurement on the sample image, select *Video* and right click in the image to activate the window context menu. Select **Imposition**. Check SHOW DATA PPOINTS box.



The profile is displayed in 2D in *SpIm* such that the color scale corresponds to the Raman intensity. The curve shape of the 2D map illustrates the shift of the Silicon position along the line of analysis that is displayed on the video image.





HOR

Explore the future

HORIBA


When the profile is displayed in 1D, peak fitting can be applied to the image in order to map the peak positions. Select **Add peak**. Click at the peak position. This adds a peak to be fitted.






Click **Approx** and then **Fit**.

When finished, click **Peak**. From *Peak* window, check MAP box for P, and click **OK**.



A new window appears with the plot of the peak positions along the line. It gives an accurate idea of the strain in the silicon structure:





Spectral Profile



- Create a Raman map out of spectra that were recorded separately.
 - Open the files that are to be made into a map
 - Select the first spectrum of the new spectral profile.
 - Input a name and click New. New SpIm window with the given name is created automatically showing the selected spectrum. Corresponding Map window is also created.
 - Add other spectra one by one by selecting and clicking **Add** or **Insert** button.
 - **Add** – adds the selected spectrum to the end of the spectral file.
 - **Insert** – inserts the selected spectrum to the position selected in the Map window

Graphic Tool Panel

- Window Context Menu
- Activated by clicking right mouse button with the mouse positioned in a window.
- Graphic Tool Panel is sensitive to the characteristics of the window it is linked to. Therefore, the specific elements listed change depending on the window the menu has been called for.
- Each Graphic Tool Panel is composed of two sections, separated by a line. Elements in the top section are common to all windows. Elements in the bottom section are window specific.

Graphic Tool Panels

Format and scale	
Fast settings	
Rescale	
Swap X axis	
Image colors	
<hr/>	
Axes	
Legend	
Spectrum	
Cursor	

Window displaying individual spectra.
Spectrum
Model

Window displaying individual video images.
Video

Format and scale	
Fast settings	
Rescale	
Swap X axis	
Image colors	
<hr/>	
Axes	
Imposition	
Image	
Scale Bar	
Cursor	

Window displaying maps.
Map
Score


Format and scale	
Fast settings	
Rescale	
Swap X axis	
Image colors	
<hr/>	
Axes	
Image	
Scale Bar	
Cursor	

Graphic Tool Panels

Format and scale
Fast settings
Rescale
Swap X axis
Image colors

Axes
Legend
Spectra


Window displaying a collection of spectra in 1D.
Splm



Format and scale
Fast settings
Rescale
Swap X axis
Image colors

Axes
Image
Scale Bar


Window displaying a collection of spectra in 2D.
Splm



Format and scale
Fast settings
Rescale
Swap X axis
Image colors


Axes
Image3D

Window displaying a collection of spectra in 3D.
Splm



HORIBA ROBIN YVON


Common Elements



Format and scale

Format the window setting. Often used options are

- Mode: Select display modes out of Single or various Multi modes.
- Rel: Allow multiple spectra/images be displayed with independent scale for easy comparison.



Shortcut icons on the right tool bar
Toggle between single and selected multi display modes

Common Elements





Image colors

Select and format color palette for the window

Common Elements





Fast settings

Activate the interactive formatting from the data window itself.

Elements that can be changed are underlined.

Shortcut icon on the right tool bar




Rescale

Rescale both X and Y axes to show the full data range.
Shortcut icon on the top tool bar

Swap x axis

Reverse x-axis direction

Window Specific Elements



Axes


Format axes and grids.



Legend

- SINGLE: Legend is displayed when Single display mode in **Formatting and scale** is selected.
- MULTI: Legend is displayed when one of Multi display mode in **Formatting and scale** is selected.


Window Specific Elements



Spectrum

Select to display spectrum as a set of Lines, Points or Bars.

Format line style, color and thickness.




Cursor

Select pointer types:


- Line: single vertical line
- Cross: vertical and horizontal lines
- Level: vertical and horizontal lines with the horizontal line position locked to the intensity
- Double: two vertical lines
- Peak: three vertical lines with positions of two outer lines locked to the approximated FWHM.

Window Specific Elements




Scale bar

Format color and font of scale bar.



Image

Select to display image with or without interpolated rendering, and change the default image color.




Imposition

Overlay elements of *Maps* or *Score* windows on top of *Video* image. Check boxes of desired elements and then select the *Map* or *Score* window to overlay.

HORIBA JOURNAL TV ON

Window Specific Elements



Spectra

Format 1D rendering of Splm window




Image 3D

Format 3D rendering of Splm window

Visual Basic Script (VBS)

A scripting language to automate LabSpec

What is VBS ?


- What is VBS?
 - In order to simplify complex or repetitive tasks, LabSpec 5 can be automated, using VBS, a well known Scripting Language.
 - Visual Basic Script is an interpreted, scripting language.
 - It allows very fast implementation, high level programming.
 - It does not requires an external development environment.
- What can be automated ?
 - Motor movements
 - Data Acquisition
 - Data Treatment
 - 2D Map or 1D Profile Generation
 - Instrument calibration
 - File import/Export
 - Video

What is Different about VBS ?

- There are basically two kinds of programming languages :
 - Compiled languages (C, C++, Pascal, Delphi, etc..)
 - They are low level programming languages that require a program called compiler to convert source code to the machine language.
 - Advantages: Very powerful, flexible and fast to execute.
 - Disadvantages: Hard to code, require a compiler, not intuitive.
 - - Interpreted languages (VBS, JavaScript, Python, PHP, ASP, etc..)
 - They are high level programming languages that do NOT require a compiler.
 - A program called interpreter converts at execution time the commands to the machine language.
 - Advantages: Very fast to implement, easy to code, very intuitive
 - Disadvantages: Much more limited and slower than compiled languages, no visual interface.

How Does it Work ?

- You said VBS is slow, but I want fast measurements !
 - LabSpec automation allows VBS (easy to code, but slow) to call LabSpec C++ (hard to code, but very fast) functions !
 - This way, you can have both an easy to code scripting language and very fast functions for acquisition and treatment.



What do I Need ?

- LabSpec 5 has an integrated script editor, with auto-completion and syntax colouring.
- It includes Full pdf documentation.
- A 'Getting Started' tutorial to get familiar with VBS and LabSpec programming is also available.

The screenshot displays the LabSpec 5 software interface. The main window is titled 'LabSpec - Spline : Spline' and contains a 'Script Editor' with the following VBS code:

```
' Constant Declaration
Const MOTOR_VALUE = 0
Const MB_STATUS_BAR = 6
Const CONSTANT_SPEED = 0
Const START_VIDEO = 0
Const STOP_VIDEO = 1
Const GET_VIDEO_ID = 2

Dim Path
Dim Mode
Dim IntegrationTime
Dim AccumulationNum
Dim AcqFrom
Dim AcqTo

Mode=ACQ_SPECTRUM
IntegrationTime=1 ' 1 sec acquisition
AccumulationNum=1 ' 1 Accumulation
AcqFrom=0 ' From=To => No MultiWindows
AcqTo=0

set filesystem = createobject("Scripting.FileSystemObject")
Path = "C:\ExperimentData\"
DoTemperatureTest "Test", 30, 40, 10, 1000, 1

Private function DoTemperatureTest(ExperimentName, StartTemperat
```

Below the script editor is a 3D plot showing a spectrum with a green grid. The plot has axes for Wavelength (nm) and Intensity. At the bottom of the interface, there are control panels for 'Laser' (441.6 nm), 'Fiber' (D 3), 'Slit' (300 μm), 'Specionometer' (450 μm), and 'Spectrometer' (511.989 nm). A 'Test' button is also visible.

Overlaid on the bottom right is an 'Adobe Reader' window displaying a PDF document titled '3 - LabSpec object Documentation'. The document content includes:

- Acq**
- Description** : Start an acquisition
- Keywords** : VBS ActiveX VB integration time accumulation
- Type** : AutoVBSAct
- long Acq**(long Mode, double IntegrationTime, long AccumulationNum, dou
- Start an Acquisition
- Mode** : Acquisition Mode :
- 0 : ACQ_SPECTRUM Start a Spectrum Acquisition
- 1 : ACQ_IMAGE Start a CCD Image Acquisition
- 2 : ACQ_LABSPEC_PARAM Start Spectrum Acquisition with LabSpec Interf
- 3 : ACQ_SPECTRAL_IMAGE Start Spectral Image Acquisition with LabSpec
- 10 : ACQ_AUTO_SHOW Add to any other Acquisition Mode : Automatically :
- IntegrationTime** : Spectrum Integration time (in seconds) (ignored if ACQ_L
- AccumulationNum** : Number of spectrum accumulation (ignored if ACQ_LA

Where Can I Use ?

- Examples in Data Acquisition
 - Spectra with different lasers
 - Spectra with different gratings
 - From low power to high power (D2, D1, D0.6,.....)
 - Depth profiles with different confocal hole diameters
- Examples in Data Processing
 - Sequences
 - filter + baseline + peak fitting
 - peak fitting + plasma position correction
 - Create new functions
 - create special data smoothing functions

Where Can I Use (cont'd) ?

- Examples in Data Display
 - Create a profile from single spectra
 - 1 D profile with access to axis definition + possibility to have non-equidistant data
 - 2 D profiles (re-create maps from single spectra)
- Examples in Data Export
 - Export in specific format for other software

Visual Basic Script (VBS) One Practical Example in Detail

HORIBA JOBIN YVON

Explore the future

HORIBA

Step 1

- Objective: Automate a Linkam heating/cooling stage
 - This example will show you how easy and how fast VBS prototyping can be, even with a relatively complex project.
- Requirements:
 - Manage a temperature ramp (From T_1 , to T_2 , Step, Speed, Holding time and Experiment name)
 - Take a video snapshot before each spectrum (save image to file)
 - Do a Raman Acquisition (save spectrum to file)
 - Create a profile Intensity vs. Temperature (save profile to file)

Step 2

1. Create a function to manage the temperature ramp
 - Temperature ramps are defined by 6 parameters. Let's write a function that can handle these parameters.

Function DoTemperatureTest (ExperimentName, StartTemperature, StopTemperature, Speed, HoldingTime, TemperatureStep)

End Function

2. Setup the Linkam stage with the correct speed and holding time
 - LabSpec specific functions are identified by "LabSpec." prefix. These functions are fully documented, specific parameters are explained and a short example is provided.

'Set Heating/Cooling Speed

LabSpec.ManageTemperature CONSTANT_SPEED, Speed, 0, Speed, 0, HoldingTime

ManageTemperature

- Overview
 - Description: Manage Linkam stage parameters
 - Keywords: Cooling, heating, speed
 - Type: AutoVBSAct
- Syntax
 - **long** ManageTemperature (**long** Mode, **double** HeatingSpeed, **double** HeatingTime, **double** CoolingSpeed, **double** CoolingTime, **double** HoldingTime)
- Parameters
 - Mode: Heating and Cooling mode
 - CONSTANT_SPEED: Use constant speed (°C/min) during heating/cooling
 - CONSTANT_TIME: Use constant time (sec) for heating/cooling
 - FREE_TEMPERATURE: Free the cooling stage (other parameters are ignored)

Step 3

3. For each point, heat or cool the Linkam stage
 - The Linkam stage can be used as any other motor (spectro, hole, slit, etc.)
 - The unique command to drive all motors is “LabSpec.MoveMotor”
 - MotorName is “Temperature” and MotorPosition (°C) for point i can be calculated by: $\text{Temp}(i) = \text{StartTemperature} + i \times \text{TemperatureStep}$

‘Calculate the number of spectra

$\text{NbSpectra} = (\text{StopTemperature} - \text{StartTemperature}) / \text{TemperatureStep} + 1$

‘Heat/Cool to the first/next step

For i = 0 **To** NbSpectra – 1

 LabSpec.MoveMotor “Temperature”, StartTemperature + I * TemperatureStep,
 “”, MOTOR_VALUE

Next

MoveMotor

- Overview
 - Description: Moving a specific motor
 - Keywords: Script, step value, string index
 - TypeL AutoVBSAct
- Syntax
 - **long** MoveMotor (**LPCTSTR** *MotorName*, **double** *PositionValue*, **LPCTSTR** *PositionName*, **long** *Mode*)
- Parameters
 - *MotorName*: Motor name (see Motor Name List)
 - *PositionValue*: Value to reach
 - *PositionName*: Position name (only for named position motors such as microscope)

Step 4

4. Capture a video snapshot

- Waiting for the video to be ready may take some time, for example to move the beam splitters in ARAMIS.
- The video image is being saved to a native “*.ngv” file here, but it can also be saved as a TIFF or JPEG image

‘Take Video Image

LabSpec.Video START_VIDEO **‘Start the live video**

Do **‘Wait till video is ready**

 VideoID = LabSpec.Video (GET_VIDEO_ID)

Loop Until VideoID > 0

‘Save the live video image

LabSpec.Save VideoID, Path & ExperimentName & “\Video_” & i+1 & “_” &
CurrentTemperature & “deg.ngv”, “ngv”

LabSpec.Video STOP_VIDEO **‘Stop the live video image**

Video

- Overview
 - Description: Display video image
 - Keywords: Camera
 - Type: AutoVBSAct
- Syntax
 - **long** Video (*long Mode*)
- Parameters
 - *Mode*: Start/Stop video
 - START_VIDEO: 0, Start video
 - STOP_VIDEO: 1, Stop video

Step 5

5. Spectrum Acquisition

- Raman Acquisition is totally similar to the video acquisition

'Take Spectrum

'Start acquisition

LabSpec.Acq Mode, IntegrationTime, AccumulationNum, AcqFrom, AcqTo

'Wait till spectrum is ready

Do

SpectrumID = LabSpec.GetAcqID()

Loop Until SpectrumID > 0

'Save the spectrum

LabSpec.Save SpctrumID, Path & ExperimentName & "\Spectrum_" &
CurrentTemperature & "deg.ngs", "ngs"

Acq

- Overview
 - Description: Start an acquisition
 - Keywords: VBS ActiveX, VB, integration time, Accumulation
 - Type: AutoVBSAct
- Syntax
 - **long** Acq (**long** Mode, **double** IntegrationTime, **long** AccumulationNum, **double** From, **double** To)
- Parameters
 - Mode: Acquisition modes
 - ACQ_SPECTRUM: 0, Start a spectrum acquisition
 - ACQ_IMAGE: 1, Start a CCD image acquisition
 - ACQ_LABSPEC_PARAM: 2, Start spectrum acquisition with LabSpec interface parameters
 - ACQ_SPECTRAL_IMAGE: 3, Start spectral image acquisition with LabSpec interface parameters
 - ACQ_AUTO_SHOW: 10, Add to any other acquisition mode to automatically show acquired data
 - IntegrationTime: Spectrum integration time (sec). Ignored if ACQ_LABSPEC_PARAM mode is selected.

Step 6

6. Add each spectrum to a profile

'If the profile does not exist and the current spectrum is the first spectrum of the profile, the profile is created.

If $i = 0$ then ProfileID = LabSpec.Profile (0, 0, SpectrumID, CurrentTemperature, "Degree", "Temperature")


'If the profile exists, the current spectrum is added the profile, and the profile is updated.

If $i > 0$ then ProfileID = LabSpec.Profile (1, ProfileID, SpectrumID, CurrentTemperature, "", "")

Profile

- Overview
 - Description: Create spectral image from spectra
 - Keywords
 - Type: AutoVBSAct
- Syntax
 - **long** Profile (**long** Mode, **long** profileID, **long** SpectrumID, **double** Value, **LPCTSTR** Unit, **LPCTSTR** Label)
- Parameters
 - *Mode*: Profile modes
 - CREATE_PROFILE: 0, Create a profile with one spectrum (ProfileID is ignored)
 - ADD_TO_PROFILE: 1, Add a spectrum to the profile
 - ProfileID: ID of a previously created profile. Relevant only for ADD_TO_PROFILE mode
 - *SpectrumID*: ID of a spectrum to add to the profile
 - *Value*: Extra dimension value for the current spectrum

Results



HORIBA JOBIN YVON

Explore the future

HORIBA