



Chemistry 4631

Instrumental Analysis

Lecture 19

Electroanalytical Chemistry

Sign Convention for Electrode Potentials (IUPAC)

Electrode potential is for half-reactions written as reductions.

Determined by the actual sign obtained when coupled with SHE in a galvanic cell.

Electroanalytical Chemistry

Sign Convention for Electrode Potentials (IUPAC)

Sign of the electrode potential, E^0 ,

- is positive when the half-cell behaves spontaneously as the cathode.
- is negative when the half-cell behaves as an anode.
- is a measure of driving force for the half-reaction.

Positive sign - Cathodic (red) reaction is spontaneous.

Electroanalytical Chemistry

Relationship between Concentration and Potential



$$K = (a_C^c)(a_D^d)/(a_A^a)(a_B^b)$$

$$\Delta G = -n F E_{\text{cell}}$$

n -- # of moles of electrons in half-reaction

F -- 96,485C

Electroanalytical Chemistry

Relationship between Concentration and Potential



$$E_{\text{cell}} = E^{\circ} - \frac{RT}{nF} \ln \frac{[C]^c [D]^d}{[A]^a [B]^b}$$

R -- 8.314 J/Kmol

T -- Kelvin

n -- # of moles of electrons in half-reaction

F -- 96,485C

ln -- 2.303log

Electroanalytical Chemistry

Relationship between Concentration and Potential

At 25°C, Nernst Equation:

$$E = E^0 - \frac{0.0592}{n} \log \frac{[C]^c [D]^d}{[A]^a [B]^b}$$

Electroanalytical Chemistry

Relationship between Concentration and Potential



$$E^\circ = RT/nF \ln K$$

Notice $E_{\text{cell}} = E^\circ$ when the activities are at unity.

Using molar concentrations instead of activities is valid only in dilute solutions.

Electroanalytical Chemistry

Examples



$$E = E^0 - \frac{0.0592}{2} \log \frac{1}{[\text{Zn}^{+2}]}$$



$$E = E^0 - \frac{0.0592}{5} \log \frac{[\text{Mn}^{+2}]}{[\text{MnO}_4^-][\text{H}^+]^8}$$

Electroanalytical Chemistry

Standard Electrode Potential, E^0

is referenced to Standard Hydrogen Electrode

-- refers to the reduction reaction (equation written as reduction)



-- independent of # of moles of reactant or product



-- positive value means reaction is spontaneous with respect to hydrogen electrode.

Electroanalytical Chemistry

TABLE 22-1 Standard Electrode Potentials*

Reaction	E^0 at 25°C, V
$\text{Cl}_2(g) + 2e^- \rightleftharpoons 2\text{Cl}^-$	+1.359
$\text{O}_2(g) + 4\text{H}^+ + 4e^- \rightleftharpoons 2\text{H}_2\text{O}$	+1.229
$\text{Br}_2(aq) + 2e^- \rightleftharpoons 2\text{Br}^-$	+1.087
$\text{Br}_2(l) + 2e^- \rightleftharpoons 2\text{Br}^-$	+1.065
$\text{Ag}^+ + e^- \rightleftharpoons \text{Ag}(s)$	+0.799
$\text{Fe}^{3+} + e^- \rightleftharpoons \text{Fe}^{2+}$	+0.771
$\text{I}_3^- + 2e^- \rightleftharpoons 3\text{I}^-$	+0.536
$\text{Cu}^{2+} + 2e^- \rightleftharpoons \text{Cu}(s)$	+0.337
$\text{Hg}_2\text{Cl}_2(s) + 2e^- \rightleftharpoons 2\text{Hg}(l) + 2\text{Cl}^-$	+0.268
$\text{AgCl}(s) + e^- \rightleftharpoons \text{Ag}(s) + \text{Cl}^-$	+0.222
$\text{Ag}(\text{S}_2\text{O}_3)_2^{3-} + e^- \rightleftharpoons \text{Ag}(s) + 2\text{S}_2\text{O}_3^{2-}$	+0.010
$2\text{H}^+ + 2e^- \rightleftharpoons \text{H}_2(g)$	0.000
$\text{AgI}(s) + e^- \rightleftharpoons \text{Ag}(s) + \text{I}^-$	-0.151
$\text{PbSO}_4(s) + 2e^- \rightleftharpoons \text{Pb}(s) + \text{SO}_4^{2-}$	-0.350
$\text{Cd}^{2+} + 2e^- \rightleftharpoons \text{Cd}(s)$	-0.403
$\text{Zn}^{2+} + 2e^- \rightleftharpoons \text{Zn}(s)$	-0.763

*See Appendix 3 for a more extensive list.

Electroanalytical Chemistry

Standard Electrode Potential, E^0

Limitations

- Use concentrations not activities (ionic strength not taken into account).
- Is temperature dependent, must specify.
- Other side reactions
- Formal potential, $E^{0'}$ - ratio of reactants and products must be unity and other[]'s are exactly known.

Electroanalytical Chemistry

Ohmic Potential or IR drop

Cells resist the flow of charge.

Ohm's Law $V = IR$ or $E = IR$

- I -- current in amper(A),
- R -- resistance in ohms(Ω)

$$E_{\text{cell}} = E_{\text{cathode}} - E_{\text{anode}} - IR$$

Electroanalytical Chemistry

Ohmic Potential or IR drop

What leads to IR drop?

- **Electrochemical cells have points of resistance throughout the cell, i.e. interfaces, wires, junctions, solution, etc...**
- **A higher potential than the thermodynamic potential, must then be applied in order to generate a current in the cell.**

Electroanalytical Chemistry

Polarization

Many electrochemical methods are based on current-voltage curves.

Use:

- Ideal polarized electrodes
- Ideal nonpolarized electrodes

Can study an electrode by coupling it with an electrode that does not readily polarize (large area, rapid reaction, reversible).

Electroanalytical Chemistry

Polarization




Figure 22-6 Current-voltage curves for an ideal (a) polarized and (b) nonpolarized electrode. Dashed lines show departure from ideal behavior by real electrodes.

Electroanalytical Chemistry

Polarization




Figure 22-7 Current-voltage curve for a cell showing ideal nonpolarized behavior between A and B (solid line) and polarized behavior (dashed line).

Electroanalytical Chemistry

Polarization Effects

- From Ohm's Law if i is plotted versus E_{app} the plot should yield a straight line with a Slope = $-1/R$.
- However, the plot for a cell deviates from a straight line at higher currents. This nonlinear behavior for a cell is called polarization. (Π)
- So a higher than expected potential must be applied to get a certain current.
- Degree of polarization = overpotential

$$\eta = E - E_{\text{eq}}$$

(E - actual potential, E_{eq} – thermodynamic potential)

Electroanalytical Chemistry

Polarization Effects

Two types of polarization

- Concentration
- Kinetic

Electroanalytical Chemistry

Polarization Effects


Concentration Polarization

Electron transfer occurs at the interface of an electrode, leaving only a thin film of solution in contact with the electrode.

- Thickness of this solution is only a few nm's and contains a finite number of molecules.
- These molecules must be replaced during the reaction.
- If the transfer of these molecules is slow, then a concentration polarization occurs.
 - How are these molecules replaced?
 - Diffusion, Migration and Convection.

Electroanalytical Chemistry

Polarization Effects



Electroanalytical Chemistry

Mass Transfer

For faradaic current to continue the species must be continuously transferred from the bulk of the solution to the electrode surface.

This occurs by:

- Convection – mechanical motion of the solution.**
- Migration – movement of ions by electrostatic attraction.**
- Diffusion – motion of species due to a concentration gradient.**

Electroanalytical Chemistry

Diffusion

- Molecules or ions move from a more concentrated region to more dilute region to eliminate any concentration gradient.
- Rate of diffusion is proportional to the concentration difference.

$$E_{cathode_applied} = E^0_{X^{+n}} - \frac{0.0592}{n} \log \frac{1}{[X^{+n}]_0}$$

- $[X^{+n}]_0$ -- conc at the electrode
- As $E_{cathode}$ becomes more negative $[X^{+n}]_0 \downarrow$ and rate of diffusion \uparrow and current \uparrow .

Electroanalytical Chemistry

Migration

Movement of ions under an electric field.

Negative ions migrate toward positive electrode and positive ions migrate toward negative electrode.

Rate of migration \uparrow as electrode $E \uparrow$.

Electroanalytical Chemistry

Convection

Mechanical transport or movement of ions or molecules in solution.

Forced convection; stirring

- decrease thickness of diffusion layer and ↓ concentration polarization.**

Electroanalytical Chemistry




Polarization Effects





Kinetic Polarization

Measure of the Rate of Electron Transfer

- Magnitude of current is limited by the rate of one or both of the electrode reactions.
- Decreased by
 - evolved gases
 - polymeric films on the electrode
- Increased by
 - catalyst

- $2\text{H}_2\text{O} + 2\text{e}^- \rightleftharpoons \text{H}_2 + 2\text{OH}^-$
-
- $2\text{HPO}_4^{2-} + 2\text{e}^- \rightleftharpoons \text{H}_2 + 2\text{PO}_4^{3-}$
-
- $\text{HPO}_4^{2-} + \text{OH}^- \rightleftharpoons \text{H}_2\text{O} + \text{PO}_4^{3-}$





Electroanalytical Chemistry

Types of Electroanalytical Methods




Figure 22-9 Summary of common electroanalytical methods. Quantity measured given in parentheses. (I = current, E = potential, R = resistance, G = conductance, Q = quantity of charge, t = time, vol = volume of a standard solution, wt = weight of an electrodeposited species.)

Assignment

- Read Chapter 22
- HW11 Chapter 22: 1, 5, 7, 9, and 11
- HW11 Due 3/22/24

