




Chemistry 4631

Instrumental Analysis

Lecture 25

Gas Chromatography

Instrumentation



Gas Chromatography

Instrumentation

Columns

Two types in general use:


- packed (stationary phase)
- open tubular or capillary

Columns determine selectivity and efficiency of the sample.

Gas Chromatography

Instrumentation


Columns



Gas Chromatography

Instrumentation

Columns



Gas Chromatography

Instrumentation

Column Materials

- **Column tubing supports the stationary phase.**
- **Needs to be:**
 - chemically inert (to prevent sample decomposition)
 - thermally stable
 - robust
 - flexible – form a coil to fit into oven

Gas Chromatography

Instrumentation

Column Materials

Early columns were metal (stainless steel, Cu, Al, Ni) but these were too reactive.

Glass is popular – but fragile

Fused silica columns – weak and subject to atmospheric corrosion. Used a protective outer Al sheath. These columns are used for most applications.

Gas Chromatography

Instrumentation

Columns

Column Activity leads to:

- tailing or skewed peaks**
- partial or total disappearance of a solute**


The surface properties of the column – catalytic or sorptive are important to consider to eliminate column activity.

Gas Chromatography


Instrumentation

Columns


Various surface groups attributed to the silica surface have been identified:




**Isolated
Silanol**



**Vicinal
Silanol**



**Geminal
Silanol**



Siloxane

Gas Chromatography

Instrumentation

Columns

- **Surface hydroxyl groups can act as proton donors in hydrogen bonding interactions**
 - -very strong sorptive sites for molecules with localized high electron density.
- **Surface siloxane bridge can act as**
 - -proton acceptors functioning as sorptive sites for molecules such as alcohols
- **Metal ions present on the surface of soda lime or borosilicate glasses can function as**
 - -Lewis acid sites absorbing molecules that have regions of localized high electron density such olefins, aromatic compounds, alcohols, ketones and amines. (Lewis acids – a species that accepts an electron pair)

Gas Chromatography

Instrumentation

Columns

Columns are deactivated before use. One method deactivates reactive hydroxyl groups by silanizing or silylating the inner surface of the column.

Typical derivatization reactions are:


- Silylation
- Acylation
- Alkylation
- Esterification

Gas Chromatography

Instrumentation

Columns

Silylation - It involves the replacement of an acidic hydrogen on the compound with an alkylsilyl group, for example, $-\text{SiMe}_3$. The derivatives are generally less polar, more volatile and more thermally stable.



Gas Chromatography

Instrumentation

Columns

Stationary Phases (GSC)

GSC of inorganic gases is the one area where packed columns are still used almost exclusively. Main adsorbents for packed columns are based on:

- silica**
- charcoal**
- alumina**
- molecular sieves**

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GSC)

Adsorbents for GSC

| <u>Chemical Type</u> | <u>Specific Surface Area (m²g⁻¹)</u> | <u>Pore Diameter (nm)</u> |
|----------------------|--|---------------------------|
| Silica | ~ 100 – 185 | 15 – 30 |
| Alumina | - | - |
| Carbon Black | 10 – 1000 | 1 – 2 |
| Molecular Sieves | 700 – 1000 | 0.5 – 2.0 |
| Porous Polymers | ~ 100 – 800 | 5 – 500 |

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

Stationary Phase is a liquid either on solid support (packed columns) or the column wall.

Solid Support – retains the liquid stationary phase and provides a large surface area. It must be thermally stable, mechanically strong to avoid breakage during the column packing process. Must have uniform pore and particle size.

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

Most common supports are based on diatomaceous earths which consist mainly of silica and minor metallic impurities. Other supports are PTFE beads, glass beads, charcoal, porous silica. Without a liquid phase these supports are essentially absorbants and would be the same as GSC.

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

The support is treated to reduce activity before coating by:

- acid wash, which eliminates metal impurities
- silanization with dimethyldichlorosilane (DMCS) or hexamethyldisilazane (HMDS) to eliminate surface activity of silanol groups.

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

Liquid Phase – required properties

- low vapor pressure
- thermal and chemical stability
- low viscosity
- nonreactive toward sample components
- wide operating temperatures (-80 to 450°C)
- reasonable solvent properties

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

Categories of stationary phases

- non-polar
- polar
- specialty phases

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

Hydrocarbons – non-polar – have been used as non-polar stationary phase.

Most have a high molecular weight for low volatility.

Examples: Squalane ($C_{30}H_{62}$), Apolane C87, long-chain n-alkanes, Apiezon L-not used much.

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

Alkylsilicone Phases – can be non-polar or polar. Polymers based on silicon-oxygen-silicone backbone – most widely used group of stationary phases. Differences in these stationary phases differ mainly in the degree of substitution on the silicon backbone.


Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

Standard polysiloxanes are characterized by the repeating siloxane backbone. Each silicon atom contains two functional groups. The type and amount of the groups distinguish each stationary phase and its properties.



R = CH₃ methyl
 CH₂CH₂CH₂CN cyanopropyl
 CH₂CH₂CF₃ trifluoropropyl



phenyl

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

Polyether Phases – polar phases – based on polyethylene glycol,



Marketed under the trade name, carbowax or superox.

Advantage - Average molecular weight of 20,000, so can be used at high temperatures.

Disadvantage is that it reacts with even trace amounts of oxygen at high temperatures, as well as decomposing to acetaldehyde and acetic acid.

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)

- High-temperature phases – substituted silicone-carborane copolymers.
- Stable up to 450°C - 500°C.
- Their advantage is in their low stationary bleed – for use with highly sensitive detectors.

Gas Chromatography

Instrumentation

High Temperature Phases




Fig. 3.18. A silicone-carborane copolymer used as a high temperature stationary phase.

Gas Chromatography

Instrumentation

Columns

Stationary Phases (GLC)


Chiral Phases - used to distinguish between enantiomers of optically active compounds.

α, β, γ – cyclodextrins – composed of six or more D(+)-glucose units bonded through a α -(1,4) – glycosidic linkage.

Gas Chromatography

Instrumentation

Chiral Phases



Gas Chromatography

Instrumentation

Columns

Stationary Phase Film Thickness

Stationary phase loadings were 20-30% w/w prior to 1970's.
Now loadings are in the order of 1-3%.

Film thickness, d_f , has a direct effect on the retention, sample capacity and elution temperature.

For best column efficiency, the film thickness is kept as thin as possible in order to reduce resistance to mass transfer in the stationary phase, C_s .

Gas Chromatography

Instrumentation

Columns

Column Selection

- **Non-polar phases usually exhibit superior lifetimes – so best to select least polar phase needed for a particular separation.**
- **Most difficult factor to access is the ability of a phase to produce a good separation. In theory selection is based on knowing the types of interactions.**

Gas Chromatography

Instrumentation

Columns

Column Selection

Types of interactions:

- London or dispersive forces (weak, non-specific).
Temporarily induced dipole. As molecular weight increases dispersive forces increase (electrons further from nucleus)
- Dipole-dipole interactions (polar molecules)
- Acid – Base interactions (proton sharing or electron transfer)

Gas Chromatography

Instrumentation

Columns

Column Selection

Selection process takes experience (many years of use).
Easiest to remember “like dissolves like.”

i.e.

Solvent

Alcohol (polar)

Alkane (non-polar)

Xylenes (non-polar)

Phase-Choice

Polar

Non-polar

Polar (highlights the slight difference in polarity among the isomers)

Gas Chromatography

Instrumentation

Columns

Methods to characterize stationary phases:

Rohrschneider Constants or McReynolds Constants

Use the relative retention of a solute – plot of the logarithm of the adjusted retention time versus Kovats retention index, I (linear).


Gas Chromatography

Instrumentation

Columns

Column Selection

Plot of the logarithm of the adjusted retention time versus Kovats retention index, I (linear).



Gas Chromatography

Instrumentation

Columns

Kovats assigned the retention index (I) of an n-alkane to equal a value 100 times its carbon number.

n-octane = 800

n-decane = 1000

n-dodecane = 12000, etc.

An I value of something different than n-alkane can be determined by “spiking” the mixture with n-alkanes. Plot $\log t'_r$ vs. I (or use equation) and I of solute is determined.

Gas Chromatography

Instrumentation

Columns

Equation:

$$I_x = 100Z + \frac{100(\log t'_{r,x} - \log t'_{r,z})}{\log t'_{r,z+1} - \log t'_{r,z}}$$

where,

Z - carbon # of n-alkane.

t'_{r,x} - adjusted retention time of component under consideration.

t'_{r,z} - adjusted retention time of alkane eluting before x

t'_{r,z+1} - adjusted retention time of alkane eluting after x.

Gas Chromatography

TABLE 3.4 Probes Used in the McReynolds and Rohrschneider Classifications of Liquid Phases

| Symbol | McReynolds Probe | Rohrschneider Probe | Measured Interaction |
|--------|------------------------|---------------------|--|
| X' | Benzene | Benzene | Electron density for aromatic and olefinic hydrocarbons |
| Y' | <i>n</i> -Butanol | Ethanol | Proton donor and proton acceptor capabilities (alcohols, nitriles) |
| Z' | 2-Pentanone | 2-Butanone | Proton acceptor interaction (ketones, ethers, aldehydes, esters) |
| U' | Nitropropane | Nitromethane | Dipole interactions |
| S' | Pyridine | Pyridine | Strong proton acceptor interaction |
| H' | 2-Methyl-2-pentanol | — | Substituted alcohol interaction similar to <i>n</i> -butanol |
| J' | Iodobutane | — | Polar alkane interactions |
| K' | 2-Octyne | — | Unsaturated hydrocarbon interaction similar to benzene |
| L' | 1,4-Dioxane | — | Proton acceptor interaction |
| M' | <i>cis</i> -Hydrindane | — | Dispersion–interaction |

Assignment

- Read Chapter 26
- HW14 Chapter 26: 1- 17
- HW14 Chapter 26 Due 04/03/24

- Read Chapter 27
- HW15 Chapter 27: 1-5, 11-19 & 21-23
- HW15 Chapter 27 Due 4/08/24

- Test 3 – PPT Lectures 15-21 - Friday April 5th

